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```

# Instrument Model
# Serial Number
# Ventilator Language
# Protocol Version
# Software Version 1
# Reserved Version
# Time.second
# Time.minute
# Time.hour
# Time.day
# Time.month

```

● TimeYear	year	#	# Time.year
● SystemTime		#	# -- System Time (yyyy/mm/dd 99h99m99s)
# 3.4.4 Monitored Parameters (Group Id 0x50)			
● Ppeak	cmH2O	#	# P max
● Pplateau	cmH2O	#	# P Plateau
● Pmean	cmH2O	#	# P mean
● PEEP	cmH2O	#	# PEEP/CPAP
● AutoPEEP	cmH2O	#	# AutoPEEP
● P0.1	cmH2O	#	# P0.1
● PTP	cmH2O*s	#	# PTP
● InspFlow	L/min	#	# Insp. Flow
● ExpFlow	L/min	#	# Exp Flow
● VTI	mL	#	# Insp. Volume
● VTE	mL	#	# Exp. Volume
● VTESpont	mL	#	# VT Exp spont
● ExpMinVol	L/min	#	# Vexp/min
● MVSpont	L/min	#	# MV Spont
● fTotal	b/min	#	# f total
● fSpont	b/min	#	# f spont
● TI	s	#	# Insp. time
● TE	s	#	# t Exp Pat
● IE		#	# I
● Rinsp	cmH2O/L/s	#	# R insp
● Cstat	mL/cmH2O	#	# Compliance
● RCexp	s	#	# RCexp
● RSB	1/(L*min)	#	# RSB
● VT_IBW	mL/kg	#	# VT/IBW
● VLeak%	%	#	# VLeak (%)
● Oxygen	%	#	# Oxygen
● Pinsp	cmH2O	#	# P insp
● VDaw	mL	#	# VDaw
● slopeCO2	%CO2/L	#	# slopeCO2
● Vtalv	mL	#	# Vtalv
● Valv	L/min	#	# V' alv
● VDawVTE	%	#	# VDaw/VTE
● PetCO2	mmHg	#	# PetCO2
● FetCO2	%	#	# FetCO2
● VeCO2	mL	#	# VeCO2
● ViCO2	mL	#	# ViCO2
● VCO2	mL/min	#	# V' CO2
● Pulse	1/min	# <C approved soon>	# Pulse
● SpO2	%	# <C approved soon>	# SpO2
● SpO2FiO2	%	# <C approved soon>	# SpO2/FiO2
● QISpO2	%	# <C approved soon>	# QI-SpO2
● QIPulse	%	# <C approved soon>	# QI-Pulse
● PetCO2 kPa	kPa	# <C approved soon>	# PetCO2 kPa
● BrethTime	min	# <C approved soon>	# Patient Breathing Time
● MVLeak	L/min	#	# MV Leak
● Thum	cel	# <GC approved soon>	# HPC Temperatur
● Typiece	cel	# <GC approved soon>	# HPC Temperatur ypiece
● Rate_10	1/min	#	# Rate 10min
● PetCO2_10	mmHg	#	# PetCO2 10min
● PetCO2_10 kPa	kPa	#	# PetCO2 10min
# 3.4.5 Special Monitored Parameters (Group Id 0x51)			
● VtiProx	mL	# <no display>	# Vtproxininsp
● VtiServo	mL	# <no display>	# Vti servo in insp
# 3.4.6 Controller State (Group Id 0x52)			
● SpO2emgLmt	%	# <for INTELLIVENT> <C approved soon>	# SpO2emergencyLimit
● SpO2incLmt	%	# <for INTELLIVENT> <C approved soon>	# SpO2increaseLimit
● SpO2decLmt	%	# <for INTELLIVENT> <C approved soon>	# SpO2decreaseLimit
● PEFIOSincFiO2	%	# <for INTELLIVENT> <C approved soon>	# PEFIOSincreaseFiO2Target
● PEFIOSdecFiO2	%	# <for INTELLIVENT> <C approved soon>	# PEFIOSdecreaseFiO2Target
● CO2incLmt	mmHg	# <for INTELLIVENT> <C approved soon>	# CO2increaseLimit
● CO2decLmt	mmHg	# <for INTELLIVENT> <C approved soon>	# CO2decreaseLimit
● fSpontInclmt	b/min	# <for INTELLIVENT> <C approved soon>	# fSpontincreaseLimit
● fSpontDecLmt	b/min	# <for INTELLIVENT> <C approved soon>	# fSpontdecreaseLimit
● MV%min	%	# <for INTELLIVENT> <C approved soon>	# MV%min(SafetyLimit)
● MV%max	%	# <for INTELLIVENT> <C approved soon>	# MV%max(SafetyLimit)
● Recruit		# <for INTELLIVENT> <C approved soon>	# Recruitmentactive
● Patient		# <for INTELLIVENT> <C approved soon>	# Patientisactive(PSS/activepatientcontrollerisusediffMV%)
● SBTrunning		# <for INTELLIVENT> <C approved soon>	# SBTrunning
● SpO2TgtS	%	# <for INTELLIVENT> <C approved soon>	# SpO2TargetShift
● CO2TgtS mmHg	mmHg	# <for INTELLIVENT> <C approved soon>	# CO2TargetShiftmmHg
● CO2TgtS kPa	kPa	# <for INTELLIVENT> <C approved soon>	# CO2TargetShiftkPa
● CPSS		# <for INTELLIVENT> <C approved soon>	# VentControllerPSS/activepatientcontrollerisactive
● NextOxy	%	# <for INTELLIVENT> <C approved soon>	# NextOxygen ->PlatformC
● NextOxyDelay	s	# <for INTELLIVENT> <C approved soon>	# NextOxygenDelay ->PlatformC
● NextPEEP	mmHg	# <for INTELLIVENT> <C approved soon>	# NextPEEP ->PlatformC
● NextPEEPDelay	s	# <for INTELLIVENT> <C approved soon>	# NextPEEPDelay ->PlatformC
● OxyFrozen		# <for INTELLIVENT> <C approved soon>	# OxygenFrozen ->PlatformC
● PEEPFrozen		# <for INTELLIVENT> <C approved soon>	# PEEPFrozen ->PlatformC
● PerMVFrozen		# <for INTELLIVENT> <C approved soon>	# PercentageMVFrozen ->PlatformC
● RecDuration	s	# <for INTELLIVENT> <C approved soon>	# RecruitmentDuration ->PlatformC
● RecPlanned		# <for INTELLIVENT> <C approved soon>	# RecruitmentPlanned ->PlatformC
● RecPreparing		# <for INTELLIVENT> <C approved soon>	# RecruitmentPreparing ->PlatformC
● RecRunning		# <for INTELLIVENT> <C approved soon>	# RecruitmentRunning ->PlatformC
● RecRunningH		# <for INTELLIVENT> <C approved soon>	# RecruitmentRunningHighPhase ->PlatformC
● TreatActOxy		# <for INTELLIVENT> <C approved soon>	# TreatmentActionOxygenation ->PlatformC
● TreatActVent		# <for INTELLIVENT> <C approved soon>	# TreatmentActionVentilation ->PlatformC
● CtrlPEEP_L	cmH2O	# <for INTELLIVENT> <C approved soon>	# ControllerPEEPLimitLow ->PlatformC
● CtrlPEEP_H	cmH2O	# <for INTELLIVENT> <C approved soon>	# ControllerPEEPLimitHigh ->PlatformC
● ManualStart		# <for INTELLIVENT> <C approved soon>	# ManualStartAllowed ->PlatformC
● PlanNextStart	s	# <for INTELLIVENT> <C approved soon>	# PlannedNextStart ->PlatformC
● StartCfullTm	s	# <for INTELLIVENT> <C approved soon>	# StartConditionsFulfilledTimer ->PlatformC
● StartingState		# <for INTELLIVENT> <C approved soon>	# StartingState ->PlatformC

E ratio (I E)

# 3.4.7 Special State (Group Id 0x53)			
● PVToolRunning		#	# P/VToolRunning
# 3.4.8 Active Alarms (Group Id 0x60)			
● Silence		#	# Silence
● ActiveAlm		#	# ActiveAlarms
# 以下は特別			
● ActiveAlm01		#	# ActiveAlarms01
● ActiveAlm02		#	# ActiveAlarms02
● ActiveAlm03		#	# ActiveAlarms03
● ActiveAlm04		#	# ActiveAlarms04
● ActiveAlm05		#	# ActiveAlarms05
● ActiveAlm06		#	# ActiveAlarms06
● ActiveAlm07		#	# ActiveAlarms07
● ActiveAlm08		#	# ActiveAlarms08
● ActiveAlm09		#	# ActiveAlarms09
● ActiveAlm10		#	# ActiveAlarms10
● ActiveAlm11		#	# ActiveAlarms11
● ActiveAlm12		#	# ActiveAlarms12
● ActiveAlm13		#	# ActiveAlarms13
● ActiveAlm14		#	# ActiveAlarms14
● ActiveAlm15		#	# ActiveAlarms15
● ActiveAlm16		#	# ActiveAlarms16
● ActiveAlm17		#	# ActiveAlarms17
● ActiveAlm18		#	# ActiveAlarms18
● ActiveAlm19		#	# ActiveAlarms19
● ActiveAlm20		#	# ActiveAlarms20
# 3.4.11 Control Settings (Group Id 0x70)			
● sModeID		# <G->15,16,20 C->3,4,15,17,18,19,20 unapproved	# ModeId
● sMode		#	# ModeName
● sApBackup		#	# ApneaBackupVentilation
● sStandby		#	# Standby
● sRate	b/min	#	# Rate
● sVT	mL	#	# TidalVolume
● sPauseTime	%	#	# PauseTime
● sFlowPat		#	# FlowPattern
● sPressTrig	cmH2O22	#	# PressureTrigger
● sPEEP	cmH2O	#	# PEEP/CPAP
● sPlow	cmH2O	#	# Plow
● sPSupport	cmH2O	#	# PressureSupport
● sFIO2	%	#	# Oxygen
● sASV	L/min	#	# ASVMV
● sPASV_Lmt	cmH2O	#	# Plimit(ASV)
● sVTsp	mL	#	# VTSetPoint
● sRRsp	b/min	#	# RRSetPoint
● sRRimv	b/min	#	# RRimv
● sPinsp	cmH2O	#	# Pinsp
● sP_Control	cmH2O	#	# PControl
● sPhigh	cmH2O	#	# Phigh
● sTrig		#	# Trigger
● sFlowTrig	l/min	#	# Flowtrigger
● sIE		#	# I
● sPeakFlow	L/min	#	# PeakFlow
● sPawPaux		#	# PaworPaux
● sETS	%	#	# ETS
● sRamp	Ms	#	# Ramp
● sBodyWt	Kg	#	# BodyWt(IBW)
● sMinVol	%	#	# %MinVol
● sTlow	s	#	# Tlow
● sVTkg	mL/Kg	#	# VT/kg
● sThigh	s	#	# Thigh
● sTlMax	s	#	# TlMax
● sTi	s	#	# Ti
● sTip	s	#	# Tip
● sPatGender		# <C->0 only C2/C3	# PatientGender
● sPatGroup		#	# PatientGroup
● sPatHeight	cm	#	# Patientheight
● sARDS		# <for INTELLIVENT> <C approved soon>	# ARDS
● sCOPD		# <for INTELLIVENT> <C approved soon>	# COPD
● sBrainInj		# <for INTELLIVENT> <C approved soon>	# Braininjury
● sHemo		# <for INTELLIVENT> <C approved soon>	# Hemodunstable(obsolete) ->No Platform
● sPneumo		# <for INTELLIVENT> <C approved soon>	# PneumoThorax(obsolete) ->No Platform
● sQuickWean		# <for INTELLIVENT> <C approved soon>	# QuickWean
● sGasSource		#	# GasSource
● sTRC		# <C approved soon>	# TRC
● sTubeSize	mm	# <C approved soon>	# Tubediameter
● sCompensate	%	# <C approved soon>	# Compensation
● sBaseflow	L/min	#	# Baseflow ->No Platform
● sSighAct		#	# Sighactivation
● sAutoRecMane		# <for INTELLIVENT> <C approved soon>	# AutomaticRecruitmentManeuver
● sUpPEEPLmt	cmH2O	# <for INTELLIVENT> <C approved soon>	# UpperPEEPlimit
● sLoPEEPLmt	cmH2O	# <for INTELLIVENT> <C approved soon>	# LowerPEEPlimit
● sDyshaemo		# <G approved soon>	# Dyshaemoglobin(obsolete) ->No Platform
● sSpO2Sensor		#	# SpO2SensorType
● sSpO2almDelay		# <C approved soon>	# SpO2alarmdelay
● sVentCtrlSts		# <for INTELLIVENT> <C approved soon>	# VentControllermanualstatus
● sPEEPCtrlSts		# <for INTELLIVENT> <C approved soon>	# PEEPControllermanualstatus
● sFIO2CtrlSts		# <for INTELLIVENT> <C approved soon>	# FIO2Controllermanualstatus
● sASV+aktiv		#	# ASV+aktiv(obsolete) ->No Platform
● sASV+VTIBW	mL/kg	#	# ASV+VT/IBW(obsolete) ->No Platform
● sASV+VTMax		#	# ASV+VTMaxenable(dobsolete) ->No Platform
● sASV+PsMV0	cmH2O	#	# ASV+PsupMV0(obsolete) ->No Platform
● sASV+PsMV100	cmH2O	#	# ASV+PsupMV100(obsolete) ->No Platform
● sASV+PsMV200	cmH2O	#	# ASV+PsupMV200(obsolete) ->No Platform
● sASV+PsMV350	cmH2O	#	# ASV+PsupMV350(obsolete) ->No Platform
● sASV+TiMax	s	#	# ASV+TiMax(obsolete) ->No Platform

● sApBackupSet		#	# ApneaBackupSetting
● sFlow	L/min	#	# Flow
● sSpeakValve	%	#	# SpeakingValve ->PlatformC
# 3.4.12 Alarm Lmts Settings (Group Id 0x71)			
● sPress_L	cmH2O	#	# LowPressure
● sPress_H	cmH2O	#	# HighPressure
● sMVe_L	L/min	#	# LowExpMinVol
● sMVe_H	L/min	#	# HighExpMinVol
● sRate_L	b/min	#	# LowRate
● sRate_H	b/min	#	# HighRate
● sVT_L	mL	#	# VTlow
● sVT_H	mL	#	# VThigh
● sLeak	%	#	# Leak
● sApneaTime	s	#	# Apneatime
● sPetCO2_L	mmHg	#	# LowPetCO2
● sPetCO2_H	mmHg	#	# HighPetCO2
● sSpO2_L	%	# <C approved soon>	# LowSpO2
● sSpO2_H	%	# <C approved soon>	# HighSpO2
● sOxy_L	%	#	# LowOxygen
● sOxy_H	%	#	# HighOxygen
● sPulse_L	1/min	# <C approved soon>	# LowPulse
● sPulse_H	1/min	# <C approved soon>	# HighPulse
# 3.4.13 GUI Unit Settings (Group Id 0x72)			
● sUnitCO2press		#	# UnitCO2pressure
● sUnitLength		#	# Unitlength
# 3.3.14 Quick Wean Settings (Group Id 0x73)			
● qQuickWean		# <for INTELLIVENT> <C approved soon>	# QuickWean
● qSaPEEP	cmH2O	# <for INTELLIVENT> <C approved soon>	# TostartSBTAdult-PEEP
● qSaOxy	%	# <for INTELLIVENT> <C approved soon>	# TostartSBTAdult-Oxygen
● qSaVTEIBW	mL/kg	# <for INTELLIVENT> <C approved soon>	# TostartSBTAdult-VTE/IBW
● qSaDelay	s	# <for INTELLIVENT> <C approved soon>	# TostartSBTAdult-Delaytime/
● qSaRSBi	1/(L*min)	# <for INTELLIVENT> <C approved soon>	# TostartSBTAdult-RSBI
● q aPEEP	cmH2O	# <for INTELLIVENT> <C approved soon>	# SBTsettingsAdult-PEEP
● q aMV	%	# <for INTELLIVENT> <C approved soon>	# SBTsettingsAdult-%MinVol
● qAaRate	%	# <for INTELLIVENT> <C approved soon>	# ToabortSBTAdult-Rateinc
● qAaOxy	%	# <for INTELLIVENT> <C approved soon>	# ToabortSBTAdult-Oxygen
● qAaPetCO2	mmHg	# <for INTELLIVENT> <C approved soon>	# ToabortSBTAdult-PetCO2inc
● qAaPetCO2_kPa	kPa	# <for INTELLIVENT> <C approved soon>	# ToabortSBTAdult-PetCO2inckPa
● qAaToITime	s	# <for INTELLIVENT> <C approved soon>	# ToabortSBTAdult-Tolerancetime
● qAaMaxDur	min	# <for INTELLIVENT> <C approved soon>	# ToabortSBTAdult-MaxDuration
● qSpPEEP	cmH2O	# <for INTELLIVENT> <C approved soon>	# TostartSBTPediatric-PEEP
● qSpOxy	%	# <for INTELLIVENT> <C approved soon>	# TostartSBTPediatric-Oxygen
● qSpVTEIBW	mL/kg	# <for INTELLIVENT> <C approved soon>	# TostartSBTPediatric-VTE/IBW
● qSpDelay	s	# <for INTELLIVENT> <C approved soon>	# TostartSBTPediatric-Delaytime/
● q pPEEP	cmH2O	# <for INTELLIVENT> <C approved soon>	# SBTsettingsPediatric-PEEP
● q pMV	%	# <for INTELLIVENT> <C approved soon>	# SBTsettingsPediatric-%MinVol
● qApRate	%	# <for INTELLIVENT> <C approved soon>	# ToabortSBTPediatric-Rateinc
● qApOxy	%	# <for INTELLIVENT> <C approved soon>	# ToabortSBTPediatric-Oxygen
● qApPetCO2	mmHg	# <for INTELLIVENT> <C approved soon>	# ToabortSBTPediatric-PetCO2inc
● qApPetCO2_kPa	kPa	# <for INTELLIVENT> <C approved soon>	# ToabortSBTPediatric-PetCO2inckPa
● qApToITime	s	# <for INTELLIVENT> <C approved soon>	# ToabortSBTPediatric-Tolerancetime
● qApMaxDur	min	# <for INTELLIVENT> <C approved soon>	# ToabortSBTPediatric-MaxDuration
● qCfBefore	min	# <for INTELLIVENT> <C approved soon>	# SBTConf-BeforestartingSBT
● qCfBetween	min	# <for INTELLIVENT> <C approved soon>	# SBTConf-Between2SBT's
● qCfRate	1/min	# <for INTELLIVENT> <C approved soon>	# SBTConf-Rate
● qCfPress	mbar	# <for INTELLIVENT> <C approved soon>	# SBTConf-Pressure
● qAfter	min	# <for INTELLIVENT> <C approved soon>	# SBTTiming-Startallowedaftertimeofday
● qBefore	min	# <for INTELLIVENT> <C approved soon>	# SBTTiming-Startallowedebeforetimeofday
● qAutoSBT		# <for INTELLIVENT> <C approved soon>	# AutoSBT ->PlatformC

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